



## SECURE RESEARCH ARCHIVE **WITH WORLD WIDE ACCESS**

The intended use of OmniStore is long-term storage of data, which in various ways have been produced in research studies from any connected modality. The data file storage structure in the archive is project based and approved users get access to the projects they are active in. The archive includes a web interface for management and export of stored data files. Data can be exported either as the original file or as a pseudonymized DICOM file.

# OmniStore – Long-term data storage with wo

## Full access wherever you are

The purpose of the system is to allow users to use a medical device on site and then be able to access generated files via a secure web access. No patient data will leave the research site and all files can be pseudonymized/identified by number for research purposes.



## The system

The system will be connected to a large medical device such as the 7 Tesla MR camera (or similar) recently installed in Lund, Sweden. Users will then be able to use the system on site and all files will be locally stored in the archive attached to the unit. The super user at the site will then allow the user to access the specific files in the project via the web back home. All the pseudonymized files will be given a unique identification number that can only be identified by the project administrator that created the pseudomimization template.

## Two factor authentication

Files can be accessed via the two factor authentication. This means that the user will login on the website and then receive an SMS with a 6 digit code used to get access. Files can be locally downloaded or exported to the DICOM server or a setup SFTP. Local files will be downloaded as zip-files, while at the device site files can be locally accessed without web access. The provided key file with the original patient data will be encrypted and the admin user have access to the password for this encrypted file.

## User levels

There are three levels of users.

**Normal user:** can only access approved projects and files defined by the project admin. DICOM files can be pseudonymized. Non-DICOM files are not pseudonymized.

**Super user:** can setup a project and add users to the project and can also allow the users to access the specific files in the project via the web back home. The super user can also setup a project and add new users. The super user will be placed at the research site and can view and change all projects.

**Project admin:** can retrieve the encryption key for identification of anonymized files. The project admin is the only person in the project that will be able to go back to the original data.

## OmniStore web login

### FEATURES

- Secure web access
- File export
- File import
- Pseudonymization of DICOM files
- Map internal drive to system for faster access of files
- Allow specific users to access specific folders

## User interface

The screenshot displays the OmniStore web interface for a project named 'TestProjekt1'. The interface includes a navigation menu at the top with options like 'Projects', 'Project info', 'Inbox', and 'Transfers'. The main content area shows a tree view of the project structure, including folders for 'John Doe', 'Misc files', and 'Dicom'. A list of files is displayed on the right, with columns for 'Name', 'Created', and 'Size'. The files are listed with their names, creation dates, and sizes. At the bottom of the interface, there are buttons for 'SFTP Export', 'DICOM Export', 'Download', and 'Upload'. The Philips logo and 'Medicinsk Bildteknik' are visible in the bottom right corner.

Name	Created	Size
✓SRe00001	2015-02-19 11:53	53.7 KB
✓SRe00002	2015-02-19 11:53	53.7 KB
✓SRe00003	2015-02-19 11:53	57.8 KB
✓SRe00004	2015-02-19 11:53	51.6 KB
✓SRe00005	2015-02-19 11:53	50.8 KB
✓SRe00006	2015-02-19 11:53	56.4 KB
✓SRe00007	2015-02-19 11:53	56.4 KB
✓SRe00008	2015-02-19 11:53	50.8 KB
✓SRe00009	2015-02-19 11:53	51.7 KB
✓SRe00010	2015-02-19 11:53	50.8 KB
✓SRe00011	2015-02-19 11:53	50.8 KB
✓SRe00012	2015-02-19 11:53	50.9 KB
✓SRe00013	2015-02-19 11:53	50.2 KB
✓SRe00014	2015-02-19 11:53	50.2 KB
✓SRe00015	2015-02-19 11:53	50.2 KB
✓SRe00016	2015-02-19 11:53	50.2 KB
✓SRe00017	2015-02-19 11:53	50.2 KB
✓SRe00018	2015-02-19 11:53	50.2 KB
✓SRe00019	2015-02-19 11:53	50.2 KB
✓SRe00020	2015-02-19 11:53	50.2 KB
✓SRe00021	2015-02-19 11:53	50.2 KB
✓SRe00022	2015-02-19 11:53	50.2 KB
✓SRe00023	2015-02-19 11:53	50.2 KB



**OmniStore**  
(Main storage)

## TECHNICAL SPECIFICATIONS

### COMPATIBLE BROWSERS

IE 8+  
Firefox 38+  
Chrome 43+  
Safari 8+

### SYSTEM REQUIREMENTS

Windows: 7-Zip, WinZip, WinRAR  
Mac OS X: iZip, Stuffit Expander, BetterZip, WinZip

### EXPORT FORMATS/TYPES

SFTP Export  
DICOM Export

Note: OmniStore is intended for research and therefore is not CE-marked.



medical imaging  
technologies

Head Office  
Medical Imaging Technologies  
Porfyrvägen 5, 224 78 Lund, Sweden

[info@medimtech.org](mailto:info@medimtech.org)  
[medicalimagingtechnologies.com](http://medicalimagingtechnologies.com)